

Anti-Leukocyte Common Antigen antibody [BRA-55] (STJ16100941)

STJ16100941

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Leukocyte Common antigen is suitable for use in Immunohistochemistry, ELISA, Immunofluorescence and Western Blot research applications.
Applications	IHC-F, IHC-P, ELISA, IF, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRA-55
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P (1:10-1:20) ELISA IF WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
Uniprot ID	PTPRC_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	