

Anti-Collagenase G antibody [1D2] (STJ16100946)

STJ16100946

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Collagenase G is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA (1:500-1:5000) WB (1:2000-1:5000)
Formulation	PBS, pH 7.2
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	