

Anti-EBFP antibody [3A6] (STJ16100947)

STJ16100947

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-EBFP is suitable for use in ELISA, Western Blot and Immunoprecipitation research applications.
Applications	ELISA, WB, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3A6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	