

Anti-PTPRR protein isoforms antibody [6A6] (STJ16100949)

STJ16100949

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PTPRR protein isoforms is suitable for use in Immunohistochemistry, Western Blot and Immunoprecipitation research applications.
Applications	IHC, WB, IP
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6A6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG2
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	5801
Gene Symbol	PTPRR
Uniprot ID	PTPRR_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	