

Mouse Anti-Chicken IgG chain antibody {Unconjugated} (STJ16100953)

STJ16100953

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated mouse monoclonal anti-Chicken IgG chain secondary antibody. For use in most research applications.
Applications	ELISA, IB
Host/Source	Mouse
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	08C
Concentration	n/a
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA IB
Formulation	Lyophilized in a 10 mM ammonium
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Polytclonal IgY (Promega)
Immunogen Region	
Specificity	
Immunogen Sequence	