

## Anti-MRP9 antibody [M9I-38] (STJ16100964)

STJ16100964

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-MRP9 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
<b>Applications</b>	ICC, IHC-F, WB, FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M9I-38
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ICC IHC-F (1:20) WB (1:50) FC
<b>Formulation</b>	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">94160</a>
<b>Gene Symbol</b>	<a href="#">ABCC12</a>
<b>Uniprot ID</b>	<a href="#">MRP9_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	