

Anti-MRP9 antibody [M9II-3] (STJ16100965)

STJ16100965

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP9 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC, IHC-F, WB, FC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M9II-3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F (1:20) WB (1:50) FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	94160
Gene Symbol	ABCC12
Uniprot ID	MRP9_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	