

Anti-Pan-T-Cell antibody [CT-6] (STJ16100972)

STJ16100972

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Pan-T-Cell is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human, Chinese, Hamster

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CT-6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Serum free tissue culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	