

Anti-Fibronectin antibody [616] (STJ16100987)

STJ16100987

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fibronectin is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC, IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Human, Dog, Ferret, Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	616
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	2335
Gene Symbol	FN1
Uniprot ID	FINC_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

