

Anti-Lymphocytic Marker antibody [IPO-74] (STJ16101001)

STJ16101001

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lymphocytic Marker is suitable for use in Flow Cytometry, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	FC, IHC-F, IHC-P, IP, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IPO-74
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

