

Anti-Tenascin-C antibody [EBS-O-166] (STJ16101005)

STJ16101005

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tenascin-C is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-F, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-O-166
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3371
Gene Symbol	TNC
Uniprot ID	TENA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	

