

Anti-UACA antibody [AE-5] (STJ16101008)

STJ16101008

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-UACA is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications. |
| Applications | FC, IHC-F, IF, WB |
| Host/Source | Mouse |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | AE-5 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1k |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 55075 |
| Gene Symbol | UACA |
| Uniprot ID | UACA_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |

