

Anti-Vascular Endothelial Growth Factor antibody [VEGF-21] (STJ16101009)

STJ16101009

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vascular Endothelial Growth Factor is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IHC-P, IF
Host/Source	Mouse
Reactivity	Human, Dog, Mouse, Rabbit, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	VEGF-21
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7422
Gene Symbol	VEGFA
Uniprot ID	VEGFA_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

