

## Anti-IgA antibody [EBS-huA] (STJ16101010)

STJ16101010

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-IgA is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, FC, IHC-F, IHC-P, IP, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	EBS-huA
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">973</a>
<b>Gene Symbol</b>	<a href="#">CD79A</a>
<b>Uniprot ID</b>	<a href="#">CD79A_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	

