

Anti-Kappa Light Chain antibody [EBS-huKappa] (STJ16101014)

STJ16101014

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Kappa Light Chain is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, FC, IHC-F, IHC-P, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-huKappa
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

