

Anti-Keratin 18 antibody [C-04] (STJ16101024)

STJ16101024

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 18 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC, IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Cat, Cow, Dog, Ferret, Goat, Hamster, Horse, Human, Marmoset, Mouse, Pig, Rat, Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C-04
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3875
Gene Symbol	KRT18
Uniprot ID	K1C18_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

