

Anti-MHC I ABC antibody [EBS-O-104] (STJ16101028)

STJ16101028

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MHC I ABC is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Monkey, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-O-104
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

