

Anti-MHC II DQ antibody [EBS-O-114] (STJ16101037)

STJ16101037

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MHC II DQ is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, IHC-F, IHC-P, IF, IP
Host/Source	Mouse
Reactivity	Human, Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-O-114
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

