

Anti-Borrelia burgdorferi antibody [EBS-I-014] (STJ16101045)

STJ16101045

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Borrelia burgdorferi is suitable for use in ELISA, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, IHC-F, IF
Host/Source	Mouse
Reactivity	Borrelia

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-I-014
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

