

Anti-Campylobacter jejuni antibody [EBS-I-001] (STJ16101046)

STJ16101046

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Campylobacter jejuni is suitable for use in ELISA, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, IHC-F, IF
Host/Source	Mouse
Reactivity	C.Jejuni

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-I-001
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgMk
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

