

Anti-Campylobacter jejuni antibody [EBS-I-002] (STJ16101047)

STJ16101047

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Campylobacter jejuni is suitable for use in ELISA, Immunofluorescence and Immunohistochemistry research applications.
Applications	ELISA, IF, IHC-F
Host/Source	Mouse
Reactivity	C.Jejuni

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-I-002
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgMk
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	