

Anti-HSV-1 gD antibody [EBS-I-042] (STJ16101055)

STJ16101055

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HSV-1 gD is suitable for use in ELISA, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, IHC-F, IHC-P, IF
Host/Source	Mouse
Reactivity	HSV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-I-042
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

