

Anti-BCL 10 antibody [EBS-T-001] (STJ16101062)

STJ16101062

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-BCL 10 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-T-001
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	8915
Gene Symbol	BCL10
Uniprot ID	BCL10_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081