

Anti-CRASH antibody [4D11] (STJ16101067)

STJ16101067

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CRASH is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4D11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	80150
Gene Symbol	ASRGL1
Uniprot ID	ASGL1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

