

Anti-MALT1 antibody [EBS-T-007] (STJ16101077)

STJ16101077

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MALT1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-T-007
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	10892
Gene Symbol	MALT1
Uniprot ID	MALT1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

