

Anti-MRP14/S100A9/calgranulin B antibody [47-8D3] (STJ16101078)

STJ16101078

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MRP14/S100A9/calgranulin B is suitable for use in Flow Cytometry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	FC, IF, IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Baboon, Cat, Cow, Dog, Goat, Guinea Pig, Human, Horse, Monkey, Pig, Rabbit, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	47-8D3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	6280
Gene Symbol	S100A9
Uniprot ID	S10A9_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

