

Anti-Cdk1/p34cdc2 antibody [POH-1] (STJ16101079)

STJ16101079

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cdk1/p34cdc2 is suitable for use in ELISA, Flow Cytometry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, FC, IF, IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Cow, Human, Mink, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	POH-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	983
Gene Symbol	CDK1
Uniprot ID	CDK1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081