

Anti-CD23 antibody [EBS-CD-016] (STJ16101099)

STJ16101099

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD23 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications. |
| Applications | FC, IHC-F, IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | EBS-CD-016 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1k |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| |
|--------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen |
| Region |
| Specificity |
| Immunogen |
| Sequence |