

Anti-CD33/SIGLE3 antibody [EBS-CD-022] (STJ16101103)

STJ16101103

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD33/SIGLE3 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-022
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF (1-2 Mu g/ml) FC (1-2 Mu g/million cells in 0 1 ml) IHC Frozen sections (1-2 Mu g/ml for 30 min at RT)
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	945
Gene Symbol	CD33
Uniprot ID	CD33_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	