

Anti-CD45 antibody [Bra55] (STJ16101108)

STJ16101108

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD45 is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, FC, IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Bra55
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
Uniprot ID	PTPRC_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	