

Anti-CD46 antibody [EBS-CD-028] (STJ16101109)

STJ16101109

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD46 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-028
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4179
Gene Symbol	CD46
Uniprot ID	MCP_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	