

Anti-CD5 antibody [Cris-1] (STJ16101110)

STJ16101110

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD5 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, IHC-F, IF, IP
Host/Source	Mouse
Reactivity	Human, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cris-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
Uniprot ID	CD5_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	