

Anti-NCAM antibody [EBS-CD-033] (STJ16101115)

STJ16101115

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, FC, IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-033
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence