

Anti-CD86 antibody [EBS-CD-042] (STJ16101126)

STJ16101126

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD86 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-042
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	942
Gene Symbol	CD86
Uniprot ID	CD86_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	