

Anti-CD100 antibody [EBS-CD-045] (STJ16101132)

STJ16101132

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD100 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Cynomolgus, Monkey, Human, Mouse, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-045
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	10507
Gene Symbol	SEMA4D
Uniprot ID	SEMA4D_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	