

Anti-CD326 antibody [EBS-CD-061] (STJ16101135)

STJ16101135

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD326 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications. |
| Applications | FC, IHC-F, IHC-P, IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | EBS-CD-061 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1k |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |