

Anti-Tumor Necrosis Factor Alpha antibody (STJ16101140)

STJ16101140

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tumor Necrosis Factor Alpha is suitable for use in Flow Cytometry, ELISA and Western Blot research applications.
Applications	FC, ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
Uniprot ID	TNFA_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	