

Anti-ICAM-1 antibody {Biotin} (STJ16101145)

STJ16101145

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ICAM-1 is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Baboon, Chimpanzee, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3383
Gene Symbol	ICAM1
Uniprot ID	ICAM1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	