

Anti-TNF-Alpha antibody (STJ16101171)

STJ16101171

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNF-Alpha is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IHC-P, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
Uniprot ID	TNFA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	