

Anti-Nitric oxide synthase antibody (STJ16101175)

STJ16101175

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nitric oxide synthase is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Rabbit
Reactivity	Cat, Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:1 000-1:5 000)
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	