

Goat Anti-Human IgM antibody {FITC} (STJ16101205)

STJ16101205

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated goat polyclonal anti-Human IgM secondary antibody. For use in most research applications.
Applications	IHC, IF, EM, WB, ELISA, FACS
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	FITC
Purification	
Dilution Range	IHC 1:40 IF EM WB ELISA FACS
Formulation	0.01 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	IgM
Immunogen Region	
Specificity	
Immunogen Sequence	