

Goat Anti-Human IgA antibody {FITC} (STJ16101208)

STJ16101208

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated goat polyclonal anti-Human IgA secondary antibody. For use in most research applications.
Applications	IHC, WB, ELISA, IF, FACS, EM
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	FITC
Purification	
Dilution Range	IHC 1:200 WB ELISA IF FACS EM
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	973
Gene Symbol	CD79A
Uniprot ID	CD79A_HUMAN
Immunogen	IgA
Immunogen Region	
Specificity	
Immunogen Sequence	