

Goat Anti-Human IgA antibody {HRP} (STJ16101209)

STJ16101209

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated goat polyclonal anti-Human IgA secondary antibody. For use in most research applications.
Applications	WB, ELISA, IHC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	HRP
Purification	
Dilution Range	WB ELISA 1:1000 IHC
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	973
Gene Symbol	CD79A
Uniprot ID	CD79A_HUMAN
Immunogen	IgA
Immunogen Region	
Specificity	
Immunogen Sequence	