

Goat Anti-Human IgG/IgM/IgA antibody {FITC} (STJ16101211)

STJ16101211

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated goat polyclonal anti-Human IgG/IgM/IgA secondary antibody. For use in most research applications.
Applications	WB, ELISA, IHC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	FITC
Purification	
Dilution Range	WB ELISA IHC
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	IgG, IgM and IgA
Immunogen Region	
Specificity	
Immunogen Sequence	