

## Rabbit Anti-Rat IgG/IgM/IgA antibody {FITC} (STJ16101232)

STJ16101232

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	FITC-conjugated rabbit polyclonal anti-Rat IgG/IgM/IgA secondary antibody. For use in most research applications.
<b>Applications</b>	WB, ELISA, IHC, IF, FACS, EM
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.2 mg/ml
<b>Conjugation</b>	FITC
<b>Purification</b>	
<b>Dilution Range</b>	WB ELISA IHC 1:200 IF FACS EM
<b>Formulation</b>	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	IgG, IgM, IgA
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	