

Rabbit Anti-Guinea Pig IgG/IgM/IgA antibody {FITC} (STJ16101235)

STJ16101235

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated rabbit polyclonal anti-Guinea Pig IgG/IgM/IgA secondary antibody. For use in most research applications.
Applications	WB, ELISA, IHC, IF, FACS, EM
Host/Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	FITC
Purification	
Dilution Range	WB ELISA IHC 1:40 IF FACS EM
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	IgG, IgM, IgA
Immunogen Region	
Specificity	
Immunogen Sequence	