

Rabbit Anti-Sheep IgG/IgM/IgA antibody {HRP} (STJ16101239)

STJ16101239

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated rabbit polyclonal anti-Sheep IgG/IgM/IgA secondary antibody. For use in most research applications.
Applications	WB, ELISA, IHC
Host/Source	Rabbit
Reactivity	Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.2 mg/ml
Conjugation	HRP
Purification	
Dilution Range	WB ELISA 1:1000 IHC
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	IgG, IgM, IgA
Immunogen Region	
Specificity	
Immunogen Sequence	