

## Rabbit Anti-Chicken IgY antibody {HRP} (STJ16101248)

STJ16101248

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	HRP-conjugated rabbit polyclonal anti-Chicken IgY secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA, IHC, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Chicken

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.2 mg/ml
<b>Conjugation</b>	HRP
<b>Purification</b>	
<b>Dilution Range</b>	ELISA 1:1000 IHC WB
<b>Formulation</b>	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	IgY
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	