

Anti-Oxytocin antibody (STJ16101252)

STJ16101252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Oxytocin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Rabbit
Reactivity	Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:400-1:800)
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	