

Anti-VACHT antibody (STJ16101255)

STJ16101255

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VACHT is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-F, IHC-P
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P (1:800-1:1 200)
Formulation	Lyophilized; reconstitute in 100 µl dist. water
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6572
Gene Symbol	SLC18A3
Uniprot ID	VACHT_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	